

GNG3-guanine nucleotide binding protein (G protein), gamma 3 Gene

Reference: PTXBC015563

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	guanine nucleotide-binding protein G(I)/G(S)/G(O) subunit gamma-3; NBP gamma-3; guanine nucleotide binding protein (G protein), gamma 3; guanine nucleotide-binding protein gamma-3 subunit; G protein subunit gamma 3		

For research use only.

GENE DESCRIPTION	guanine nucleotide binding protein (G protein), gamma 3
GENE ID	2785
ORIGIN SPECIES	Human
ACCESSIONS	BC015563
VECTOR	pDONR223
NCBI SYMBOL	GNG3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaagtgagaccccggtgaacagcactatgagtattggcagcagcgaagatgtggaacagcttaagattgaagccagcttgtgtcggataaaggtgtccaagcagcagcagacctgatgacttactgtgatgccacgcctgtgaggatccctcatcaccctgtgccactt
cggagaacccctccgggagaagaagtctctgtgctctcctctga