



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

CSDA-cold shock domain protein A Gene

Reference: PTXBC015564

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

CSDA; CSDA1; DBPA; ZONAB; Y-box-binding protein 3; DNA-binding protein A; ZO-1-associated nucleic acid-binding protein; cold shock domain-containing protein A; cold-shock domain containing A1; cold-shock domain protein A; single-strand DNA-binding protein

For research use only.

GENE DESCRIPTION	cold shock domain protein A
GENE ID	8531
ORIGIN SPECIES	Human
ACCESSIONS	BC015564
VECTOR	pDONR223
NCBI SYMBOL	CSDA
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr