

HOXB9-homeobox B9 Gene

Reference: PTXBC015565

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HOX-2.5; HOX2; HOX2E; homeobox protein Hox-B9; homeo box 2E; homeo box B9; homeobox protein Hox-2.5; homeobox protein Hox-2E; homeobox B9		

For research use only.

GENE DESCRIPTION	homeobox B9
GENE ID	3219
ORIGIN SPECIES	Human
ACCESSIONS	BC015565
VECTOR	pDONR223
NCBI SYMBOL	HOXB9
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtccatttctggagcgttagcagctattatgtcactcgcataagtcacgagagtgaggacgcgcctccagccaaagtttcttctgcccagtagcgcagctcgcggcagccggccacgcggagcacctggagttccctctgtcagcttcagccaaagcgcgggtgttcggc
gcctctgtggcgcgcgtgagccgcacgcgtccgggagcctgcccgtcgtctaccaccttaccacccagccagggcgtccgcgcgcgcgagagcaggtacctccgcacctggctggagccggcgcgcgcgcgaagcggcccccgggagggccagcgcgggtgaagg
cggagccgctgtgtggcgcgcctgggagctgtctaaacagggcagcccgagtagctgttggaaacttcggcggcagggagggcgtgtgtctaatcaagaccggcgtacgggacaataaaatttgcgaaggaagcagggacaaagagagggccgatcaaccaaccct
ccgccaactggctgcagctcgtcttccggaaaaagcgtgtccctacacaaataccagacgctggagctagagaaggagtttctgttcaatgtacctcaccagggaccgtaggcacgaagtggccagactcctcaatctgagtgagagacaagtcaaatctgtttcagaaccg
gcggatgaaaatgaagaaatgaataaggagcaggcgaagagtaa