

FAM165B-family with sequence similarity 165, member B Gene

Reference: PTXBC015596

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein FAM165B; FAM165B; C21orf51; SMIM11; SMIM11B; small integral membrane protein 11A; family with sequence similarity 165, member B; small integral membrane protein 11		

For research use only.

GENE DESCRIPTION	family with sequence similarity 165, member B
GENE ID	54065
ORIGIN SPECIES	Human
ACCESSIONS	BC015596
VECTOR	pDONR223
NCBI SYMBOL	FAM165B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaattggaaggttctgagcacgtgccctgctgtatatcttggcagcaaaacattattctctgcctgacatttctgggtgaaaatgtatcaagaaaaagggtggaggcaaaacaacaaactggaggctgaaaggaagaagcaatcagagaaaaagataactga