

ELSPBP1-epididymal sperm binding protein 1 Gene

Reference: PTXBC015598

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	E12; EDDM12; EL149; HE12; epididymal sperm-binding protein 1; epididymal protein 12; epididymal secretory protein 12; epididymis luminal secretory protein 149; epididymal sperm binding protein 1		

For research use only.

GENE DESCRIPTION	epididymal sperm binding protein 1
GENE ID	64100
ORIGIN SPECIES	Human
ACCESSIONS	BC015598
VECTOR	pDONR223
NCBI SYMBOL	ELSPBP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacccgatgggtccagttacctgttgggatggacaaccttcttctattctatgagtcagtgaggatgcatgaggatgtgtcttcttccactacaaggatctgttacttcaactgcacccatattcatagcttatcccttgggtgtgccaccagagccgtgtacaacggccagtgg
aagtactgccagagtgaaattaccacacgtgtatcttcccttcatctatcgaggaaagcttataacagctgcatctcccagggcagcttcttaggcagctctgtgtgctcagtcacctctgtcttcgatgagaacagcagtggaattctgtgaaacgaatgagtatggggaaattcttc
agggaagcctgcatcttccctccatctacagaataatgtgtctctgattgcatggagatgaaagcaacaagctctgtgtcccaaccagagaacatggataagatggaaagtggagttctgtgccgacaccagaattccgctgtgtccctggccttctgtcactttccgttcaa
ctataaaacaagaattattttaactgcactaacaaggatcaaaaggagaacctgtgtgtgtgcaacttctacaactacgaccaagaccacacctgggtgtattgctga