



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

MAPK1IP1L-mitogen-activated protein kinase 1 interacting protein 1-like Gene

Reference: PTXBC015621

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C14orf32; MISS; c14_5346; MAPK-interacting and spindle-stabilizing protein-like; MAPK-interacting and spindle-stabilizing protein; mitogen-activated protein kinase 1 interacting protein 1 like		

For research use only.

GENE DESCRIPTION	mitogen-activated protein kinase 1 interacting protein 1-like
GENE ID	93487
ORIGIN SPECIES	Human
ACCESSIONS	BC015621
VECTOR	pDONR223
NCBI SYMBOL	MAPK1IP1L
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctgatgaattctgttgcagatgcactacctgaacacccccctgcaaaacctctctgtgagcaatacaaaaacctggccaacctctcaaggctggccaggctccaacctctggaataatccgagtgctccaattcagtgccattctggactccacaaagtgcaacacctcactgt
gctctttggaccagcaaccaaggaaatgtatccctccgtgctctccaccggaccacctccaggaccctccagacccttctctctccggaccatcatgtccccacctgtgtgtgctctatccagccccactgtgctggcccttggccccacaggggccatatctcaacaaatgatgccctt
tccagagctaccacagaccatagtgtgctaccacagatccagctgcagctggctctttagcttcatggggatccatgtctcttggaaccttggggccaggaaatggaggggcagtatcttaccctaataatgcatcatccatccaggccatataccctctctctctccaggcccccttg
ggcagcaccacctgttcatggggcaccgttccaccaggagcctggggaccaccagaccatatctgcccctacaggaatgataccacacaggactctactctccagtaatactttcaagtgccttcaggacctctgtgtgtcccaaatgctgtgtgccccattctacat
taa

More info: Please contact us for any question: contact@proteogenix.fr