

TRIML1-tripartite motif family-like 1 Gene

Reference: PTXBC015684

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RNF209; RING finger protein 209; tripartite motif family like 1		
For research use only.			

GENE DESCRIPTION	tripartite motif family-like 1
GENE ID	339976
ORIGIN SPECIES	Human
ACCESSIONS	BC015684
VECTOR	pDONR223
NCBI SYMBOL	TRIML1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaggagatgctaagaaattcagcacggagataacgctggaccagccacagctaatacctatctctgtttctcggaggatctgaagagtgtgaaatatgggggaagcagacagcagctacccgacaacccggaaagatttgaccagtctgcgactgtgctgggtactcagatcttc
accagtgggagacactactgggaggtggaggtgggaaacaagaccgagtggaagtgggcactctgcaaggactctgtgagcagaagggaatctcccaagccactggggacctgttctactaatagggttaaaatcggagatgattacagcctctgggtctcgtcactttgaaa
ggtcagcacgtcagagagcctgtgtgtaagggtgtgtctctcctggactatgaatctggacatatacattctacaacgggacggatgaatccctcatctacagcttcccagggctcttccaagagccctcaggcctatctttcccccctgcctcccaatgaggggacaaacacagacc
ctctcaccatctgctcaactgaacagccacgtctga