

RAB22A-RAB22A, member RAS oncogene family Gene

Reference: PTXBC015710

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RAB22A, member RAS oncogene family; GTP-binding protein RAB22A; ras-related protein Rab-22A; rab-22		

For research use only.

GENE DESCRIPTION	RAB22A, member RAS oncogene family
GENE ID	57403
ORIGIN SPECIES	Human
ACCESSIONS	BC015710
VECTOR	pDONR223
NCBI SYMBOL	RAB22A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgctgaggagctcaaatgtgtctgctcggggatagcaggtgttaggtaaatcagatattgtgtggcggtttgtggaagacagttttgatccaacatcaaccaacaataggggcatctttatgaccaagactgtccagtacaaaatgagctacataaattcctaattctgggatacagc
tggacaagaacgatttcgtgcttagcaccatgtactatcgagggtcggctgcagctataatcgtttatgatacacaagaagagacattttcaacattaaagaattgggtgaaagagcttcgacagcatggcccacctaattgtatgtccattgcaggaaataaatgtgatcttatcgat
gtaagagaagtcattggagagatgcaaggactacgccgactctattcatgcaatttttagagaccagcgcaaaaaacgcgataaacataaatgaactctttatagaattagtcgaagaattccactgacgccaacctgccatctggcggttaagggtctcaactccgaagaca
gccttcagagccaagcgagctgctgctga