

CXCL2-chemokine (C-X-C motif) ligand 2 Gene

Reference: PTXBC015753

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CINC-2a; GRO2; GROb; MGSA-b; MIP-2a; MIP2; MIP2A; SCYB2; C-X-C motif chemokine 2; GRO2 oncogene; MGSA beta; MIP2-alpha; chemokine (C-X-C motif) ligand 2; gro-beta; growth-regulated protein beta; macrophage inflammatory protein 2-alpha; melanoma growth sti		

For research use only.

GENE DESCRIPTION	chemokine (C-X-C motif) ligand 2
GENE ID	2920
ORIGIN SPECIES	Human
ACCESSIONS	BC015753
VECTOR	pDONR223
NCBI SYMBOL	CXCL2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcccgccacgctctccgccccccagcaatccccggctctcgccggtggcgctgctctctctctggtggccgcccagccggcgccgagcaggagcgccccctggccactgaactgcgctgccagtcttgcagaccctgcaggggaattcacctcaagaacatccaaagt
gtgaagggtgaagtcctcccgacccactgcgccccaaaccgaagtcatagccacactcaagaatggcgagaaagctgtctcaacccgcacgcccatggttaagaaaatcatcgaaaagatgctgaaaaatggcaaatccaactga