



Reference: PTXBC015758

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	170A; DDX13; HLP; SKI2; SKI2W; SKIV2; SKIV2L1; THES2; helicase SKI2W; SKI2 homolog, superkiller viralicidic activity 2-like; helicase-like protein; superkiller viralicidic activity 2-like; Ski2 like RNA helicase		

For research use only.

GENE DESCRIPTION	superkiller viralicidic activity 2-like (S. cerevisiae)
GENE ID	6499
ORIGIN SPECIES	Human
ACCESSIONS	BC015758
VECTOR	pDONR223
NCBI SYMBOL	SKIV2L
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr