

SH2B3-SH2B adaptor protein 3 Gene

Reference: PTXBC015786

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	IDDM20; LNK; SH2B adapter protein 3; lymphocyte-specific adapter protein Lnk; signal transduction protein Lnk; SH2B adaptor protein 3		

For research use only.

GENE DESCRIPTION	SH2B adaptor protein 3
GENE ID	10019
ORIGIN SPECIES	Human
ACCESSIONS	BC015786
VECTOR	pDONR223
NCBI SYMBOL	SH2B3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagtttctgtgaacccggctggagtcagtggtgcaatctggctcactgcaactccacttctctggttcaagcgattctgctcgacctctcaagtagctgggattacaagcaccagccaccatgcctgctaattttgtatttttag