

RAB38-RAB38, member RAS oncogene family Gene

Reference: PTXBC015808

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RAB38, member RAS oncogene family; NY-MEL-1; rrGTPbp; ras-related protein Rab-38; Rab-related GTP-binding protein; melanoma antigen NY-MEL-1		

For research use only.

GENE DESCRIPTION	RAB38, member RAS oncogene family
GENE ID	23682
ORIGIN SPECIES	Human
ACCESSIONS	BC015808
VECTOR	pDONR223
NCBI SYMBOL	RAB38
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcaggccccgcacaaggagcacctgtacaagttgctgtgattggcgacgtggcggtgggaagaccagatcatcaagcgctacgtgcaccagaactctcctgcactaccgggccacaatcggtggaacttcgcgctcaaggtgctccactgggacccggagactgtgtgc
gcctgcagctctgggafatcgaggtcaagaaagatttggaaacatgacgaggggtctattaccgagaagctatgggtgcaattattgtcttcgatgcaccagggccagccacattgaaagcagtggaaggtggaaaaatgattggactccaagttaagtcctaatggcaaccggttc
agtgtttttgtggccaacaatgtgaccaggggaaggatgtctcatgaacaatggcctcaagatggaccagttctgcaaggagcacgggttcgtaggatggttgaacatcagcaaggaaaataaacaattgatgaagcctccagatgcctgggtgaacacatacttgaatgagtg
tgacctaatggagtctattgagccggagctgtgaagcccatctcacatcaaccaaggtggccagctgtctgctgtgccaatcctag