

MAP1LC3A-microtubule-associated protein 1 light chain 3 alpha Gene

Reference: PTXBC015810

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ATG8E; LC3A; MAP1ALC3; MAP1BLC3; microtubule-associated proteins 1A/1B light chain 3A; MAP1 light chain 3-like protein 1; MAP1A/1B light chain 3 A; MAP1A/MAP1B LC3 A; MAP1A/MAP1B light chain 3 A; autophagy-related ubiquitin-like modifier LC3 A; microtubul		

For research use only.

GENE DESCRIPTION	microtubule-associated protein 1 light chain 3 alpha
GENE ID	84557
ORIGIN SPECIES	Human
ACCESSIONS	BC015810
VECTOR	pDONR223
NCBI SYMBOL	MAP1LC3A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccctcagaccgaccttcaagcagcggcgagcttcgccgaccgctgtaaggaggtacagcagatccgcgaccagcaccagcaaaatcccggtgatcatcgagcgctacaagggtgagaagcagctgccgctctggacaagaccaagttttgggtcccgaccatgtcaaca
tgagcgagttggtcaagatcatccggcgccgctgcagctgaacccacgcaggccttctctgctggtgaaccagcacagcatggtgagtggtgtccacgccatcgccgacatctacgagcaggagagaagacgaggacggcttctctatggtctacgcctccaggaaaccttcg
gcttctga