

## LGALS8-lectin, galactoside-binding, soluble, 8 Gene

Reference: PTXBC015818

|                |                                                                                                                                                                                                                        |                     |                |
|----------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                                                    | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                                                   | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | Gal-8; PCTA-1; PCTA1; Po66-CBP; galectin-8; Po66 carbohydrate binding protein; galectin-8g; lectin, galactoside-binding, soluble, 8; po66 carbohydrate-binding protein; prostate carcinoma tumor antigen 1; galectin 8 |                     |                |

**For research use only.**

---

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | lectin, galactoside-binding, soluble, 8                                                          |
| <b>GENE ID</b>          | 3964                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC015818                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | LGALS8                                                                                           |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon)

atgtgtcttaaacacctacagaatatcatctataaccggtaacccgtatgttggcaccattcccgatcagctggatcctggaacttggatgtgatgtgtggcatgttcttagtgacgcagacagattccaggtggatctgcagaatggcagcagtgtaaacctcgagccgatgtggcc  
tttcatttcaatctctgttcaaaaggccggctgcatgtttgcaatacttggataaatgaaaaatggggacgggaagagatcacctatgacacgccttcaaaagagaaaaagcttttgagatcgtgattatggtgctaaaggacaaattccaggtggctgtaaatgaaaaacatactctgctcta  
tgccacaggtcggccagagaaaaatagacactctggcatttatggcaaatgtaataatcactcaattgggttttagcttcagctcggacttacaagataccaagcatctagctggaactgacagagataagtagagaaaaatgttccaaagtctggcagccccagcttctagtaataga  
ggaggagacatttcaaaatcgacccagaactgtctacccaagagcaaatgctgactgcaatcacacttgaacttgacacaaataaccacctatgaactatgtgtcaaaagagcctgccattcgtgcaaggttgaacacccccatgggcccctggacgaactgtcgtcttaaaaggagaa  
gtgaatgcaaatgcaaaagctttaatgttgactactagcaggaataatcaaaagatattgctctacacttgaaccacgcctgaataatgaagcatttgaagaattcttctcaggagctcctgggagagaagagagaaatattacctcttccatttagctctgggatctttgagatg  
ataatttactgtgatgttagagaattcaagggttcagtaaatggcgtacacagcctggagtacaaacacagatttaagagctcagcagatttgacacgctggaataatggagacatccactactggaagtaaggagctgtag