

MRPS35-mitochondrial ribosomal protein S35 Gene

Reference: PTXBC015862

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HDCMD11P; MDS023; MRP-S28; MRPS28; 28S ribosomal protein S35, mitochondrial; 28S ribosomal protein S28, mitochondrial; MRP-S35; S28mt; S35mt; mitochondrial ribosomal protein S28; mitochondrial ribosomal protein S35		

For research use only.

GENE DESCRIPTION	mitochondrial ribosomal protein S35
GENE ID	60488
ORIGIN SPECIES	Human
ACCESSIONS	BC015862
VECTOR	pDONR223
NCBI SYMBOL	MRPS35
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcaaggagggaatctagaacttttaagattccaattttctgcaatttgactcctgttagcaattaaaaagcactgtgaagccctaaagattttgcactgagtgccagccgactggagcagtgacgagaaatgtgagaagcattttccaattgaaattgacagcactgattatgtttcatc
aggacatctgttcggaacccagagcagcagtagtagtcttaagagtaaaagctttccagtttgatttagatgacacgcaagaagaataaataaactgttaggagagcgatactgcaagaccacagatgtgcttaccatcaaacagatagggtgcccttaaggaggcagaattacgat
tatgcagtgtaactaactacagtggtataccatgagctcttgaatactgaagaatgggaaaaaagtaagactgaagcagacatggaagagtatatatgggaaatagctcatcagaagaataatcctggaacgcttccagatgaaagctgctgagaaaaatggaataataaagaa
gagctccttggtactaaagaattgaagagtacaaaagctgtgttagtcttaaaaatgaggaggaatgaaattccatttctcagtacaagaatccgtgaagagactattaaatgtgacatga