

TMEM186-transmembrane protein 186 Gene

Reference: PTXBC015912

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C16orf51; transmembrane protein 186		
For research use only.			

GENE DESCRIPTION	transmembrane protein 186
GENE ID	25880
ORIGIN SPECIES	Human
ACCESSIONS	BC015912
VECTOR	pDONR223
NCBI SYMBOL	TMEM186
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctgccctctccagctgtgctaggtttcggggaaagctgtgtgggaaagccctccatggcgtgtgtgctgcagtgggcaggaggatcctaagagtggtgggcagcagttcacccatctcgaaggagaaactacaaacgcagagactgagaaattctggatgtttacc
gttttgatgccatcagaaccttcgggttcctgtcagactgaagtggcacagactgccctgacagtggtgactttgccaccaggctattactgtactccaggccctcctcactctcaacaccgtgtgctcatgagtggtgatcgggctttgccctgacctgctgtgctggatgactattt
cttacggagactggttggtatctgtatctgaatgagctgtggcaccatgctgcgggtggccatctgaactctgggctggcggcagacacatactgtcccatggcagatgtgattccctgacagaaaccaaggaccgacctcaggagatgtttgtcgatccagcggtacagtgga
aacagacctctacgtcaccctgcgctatggacgcatcctggacagagagcgtttcacacaggtgtttgggtacatcagatctcaagtga