



Reference: PTXBC015914

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ASK2; MEKK6; mitogen-activated protein kinase kinase kinase 6; apoptosis signal-regulating kinase 2		

For research use only.

GENE DESCRIPTION	mitogen-activated protein kinase kinase kinase 6
GENE ID	9064
ORIGIN SPECIES	Human
ACCESSIONS	BC015914
VECTOR	pDONR223
NCBI SYMBOL	MAP3K6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]