



Reference: PTXBC015920

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CSNB1H; guanine nucleotide-binding protein G(I)/G(S)/G(T) subunit beta-3; G protein, beta-3 subunit; GTP-binding regulatory protein beta-3 chain; guanine nucleotide binding protein (G protein), beta polypeptide 3; guanine nucleotide-binding protein G(I)/G		

GENE DESCRIPTION	guanine nucleotide binding protein (G protein), beta polypeptide 3
GENE ID	2784
ORIGIN SPECIES	Human
ACCESSIONS	BC015920
VECTOR	pDONR223
NCBI SYMBOL	GNB3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]