



Reference: PTXBC015922

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MILT1; OIP106; trafficking kinesin-binding protein 1; 106 kDa O-GlcNAc transferase-interacting protein; O-linked N-acetylglucosamine transferase interacting protein 106; OGT(O-GlcNAc transferase)-interacting protein 106 KDa; milton homolog 1; trafficking		

GENE DESCRIPTION	trafficking protein, kinesin binding 1
GENE ID	22906
ORIGIN SPECIES	Human
ACCESSIONS	BC015922
VECTOR	pDONR223
NCBI SYMBOL	TRAK1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1