

ZNF706-zinc finger protein 706 Gene

Reference: PTXBC015925

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HSPC038; PNAS-106; PNAS-113; zinc finger protein 706		
For research use only.			

GENE DESCRIPTION	zinc finger protein 706
GENE ID	51123
ORIGIN SPECIES	Human
ACCESSIONS	BC015925
VECTOR	pDONR223
NCBI SYMBOL	ZNF706
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctcgtggacagcagaaaatcagtcagcagaaaaatgccaaaagcaagctggacaaaagaagaacaaggacatgccaaaagcgtgctgccaagctgccttaatatatctgcactgtctgttaggacacaaatgccagaccctaagaccttaagcagcactttgagagc
aagcatcctaagactccacttctccagaattagctgatgttcaggcataa