

KRCC1-lysine-rich coiled-coil 1 Gene

Reference: PTXBC015927

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CHBP2; lysine-rich coiled-coil protein 1; cryptogenic hepatitis-binding protein 2; lysine rich coiled-coil 1		

For research use only.

GENE DESCRIPTION	lysine-rich coiled-coil 1
GENE ID	51315
ORIGIN SPECIES	Human
ACCESSIONS	BC015927
VECTOR	pDONR223
NCBI SYMBOL	KRCC1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagcattcaaaagacatgatgactcttttcaagatgaactgaagattatataaagtacagaagccagaggcttagagccaaagactgtttcagaagatgaaaggggactatttgaaacctgtgggtacaaggagaggttaattccagaccacgtatagaatgttgaccaga
gactcccatctgaaacatccagacctaccaagatcatgcaatattccacaacagtggaataatcggttgctcagtggttaccagcccatgacagcagattgagactagactctctgagctactgtcagttcacgaggagactgttttcagaaaaaccagtacccctgaactttaatcaaca
agaatataattgtggctcacatggtgtagaacaatagattacaagcacttctctcagataacagtaccagttactcaagccagtcacaacagatacatcagaagaggaaaaagccagagggaaggcagagaaaaatcagaggaggagcggtctaaagcataagagaaaaaa
gctgcgaggaaattgacttagacaacaagagcatccaagaaagaaaaacagaggtggaaatagaaccgtacatgtcagtcagaaaagcttaagaatcgaaggagaaaaaagccgagatgtagtcttaagaaagaggaacgttaagcgtacaaaaagaaaaaggaacaa
ggccaagaaaggacagaggaggaatgcttgggaccagttcttcttgattttga