

RAB35-RAB35, member RAS oncogene family Gene

Reference: PTXBC015931

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RAB35, member RAS oncogene family; H-ray; RAB1C; RAY; ras-related protein Rab-35; GTP-binding protein RAY; ras-related protein rab-1c (GTP-binding protein ray)		

For research use only.

GENE DESCRIPTION	RAB35, member RAS oncogene family
GENE ID	11021
ORIGIN SPECIES	Human
ACCESSIONS	BC015931
VECTOR	pDONR223
NCBI SYMBOL	RAB35
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcccgaggactacaccactcttcaagctgctcatcggcgacagcggtgtgggcaagagcagttactgttgcgttttcagacacactttctcaggcagctacaccacgatcggagtgattcaagatccggaccgtggagatcaacggggagaaggtgaagctgcag
atctgggacacagcggggcaggagcgtctccgaccatcacctccacgtattatcgggggacccacggggcattgtgtttacgacgtcaccagtgccgagtccttgcacagtcacagcgtggcttcacgaaatcaaccagaactgtgatgtgtgcccgaatattagtggtataa
gaatgacgacctgagcggaaagtggtggagacggaagatgcctacaaatcgccggcagatgggcatccagttgttcgagaccagcgccaaggagaatgtcaacgtggaagagatgttcaactgcacacggagctggctccgagcaaaagacaacctggcaaacacg
agcagcaacaacagacgatgtgtggaagctcacgaagaacagtaaacgaaagaaacgctgctgctaa