

## RBBP9-retinoblastoma binding protein 9 Gene

Reference: PTXBC015938

|                |                                                                                                                                                                                                       |                     |                |
|----------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                                   | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                                  | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | BOG; RBBP10; B5T overexpressed gene protein; RBBP-10; RBBP-9; retinoblastoma-binding protein 10; retinoblastoma-binding protein 9; retinoma-binding protein 9; RB binding protein 9, serine hydrolase |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | retinoblastoma binding protein 9                                                                  |
| <b>GENE ID</b>          | 10741                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC015938                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | RBBP9                                                                                             |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atggctctctagcaaggcagtgattgtcccggaacggaggcggggatgtgaccacccacggctggtatggctgggtgaaaaaggagctggagaagatacctggttcagtggttgctaaaaacatgcccgaccaattacagcacgagagagcatctggctgcccttcaggag  
acagagctgcactgtgatgagaagactatcatcattggccacagttctggggccatcgccgcatgaggtatgcagaacacatcgagtatactgctattgtattagtgctgctacacatcagactgggggatgaaatgagcgtgcaagtggatactcaccgccccctggcagtgga  
gaagalcaaggccaactgcccttacctgtgagtttggtctactgacgacccgttccttcctggaaggaacaacaagaagtgccgcataggttggaaccaaattgcacaaattcactgactgtggccacttcagaacacagagttcatgaactgattactgtgttaaagctttgctgaa  
agtaccagcatag