

NIP7-nuclear import 7 homolog (S. cerevisiae) Gene

Reference: PTXBC015941

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NIP7, nucleolar pre-rRNA processing protein; 60S ribosome subunit biogenesis protein NIP7 homolog; CGI-37; HSPC031; KD93; nuclear import 7 homolog		

For research use only.

GENE DESCRIPTION	nuclear import 7 homolog (S. cerevisiae)
GENE ID	51388
ORIGIN SPECIES	Human
ACCESSIONS	BC015941
VECTOR	pDONR223
NCBI SYMBOL	NIP7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgcctttgactgaagagagacccgtgtcatgtttgagaagatagcgaatacattggggagaatctcaactgctgtggaccggccgatggcacctactgtttccgtctgcacaacgaccgggtgtactatgtgagtgagaagattatgaagctggccccaatttccgggga
caagctgtgtgctgtgggacctgtttgaaaattcactaaaaccacaagtttcggtgtgcacgtcacagctctgattaccttgaccttatgccaaagtataaagtttggataaagcctgtgtgcagagcagtccttctgtatgggaaccatgtgttgaatctgtctgtgggtcgaatcactga
aaatacttctcagtaccagggcgtgtgtgtgtactccatggcagacatcccttgggtttgggtggcagccaaatctacacaagactgcagaaaagtagaccccatggcgattgtgttattcatcaagcagacattggggaatatgtcggcatgaagagacgttgacttaa