

ORM2-orosomucoid 2 Gene

Reference: PTXBC015964

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AGP-B; AGP-B'; AGP2; alpha-1-acid glycoprotein 2; AGP 2; OMD 2; alpha-1-acid glycoprotein, type 2; orosomucoid 2		

For research use only.

GENE DESCRIPTION	orosomucoid 2
GENE ID	5005
ORIGIN SPECIES	Human
ACCESSIONS	BC015964
VECTOR	pDONR223
NCBI SYMBOL	ORM2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgctgtcctgggtttctacagctcctcctcctcctgctggaagccagatccattgtgtgccaacctagaccgggtgcccatcaccaagccacccctggaccggatcactggcaagtgggttttatcgcacggcctttcgaacgaggagtagacaataagtcggttcaggaga
tccaagcaacctcttttactttacccccacaagacagaggacacgatctttctcagagagtagacagaccgcccagaaccagtgcttctataactccagttacctgaatgtccagcgggagaaatgggaccgtctccagatacaggaggccgagaacatgttgcctcactgctgttccttag
ggacaccaagacctgatgtttggttcctacctggacgatgagaagaactggggctgtcttctatgctgacaagccagagacgaccaaggagcaactgggagagttctacgaagctctcactgctgtgtcattccagggtcagatgtcatgtacaccgactggaaaaggataagtg
agccactggagaagcagcagcagagaaggagaggaaacaggaggagggggaatcctag