

MRPS36-mitochondrial ribosomal protein S36 Gene

Reference: PTXBC015966

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DC47; MRP-S36; 28S ribosomal protein S36, mitochondrial; S36mt; mitochondrial ribosomal protein S36		

For research use only.

GENE DESCRIPTION	mitochondrial ribosomal protein S36
GENE ID	92259
ORIGIN SPECIES	Human
ACCESSIONS	BC015966
VECTOR	pDONR223
NCBI SYMBOL	MRPS36
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatgggcagcaagatggcgctgctagtaggggtcgttcaggtagtcgaaccacacactccatttaataaggtttctgacagaagagacaatcctaaaccaatgtatcagaagctttgagatcagcagggtaccatctcactcttctgtaatttcacaacattctaaaggaagtaaatcacc
agatttgcctgatgtatcaggggtccaccagacactgcagaataataaaaaacattacctcagaataacagaaggaaactgtgtctcaagaagaatggaatttatccaactggaggtcctgaataa