

S100A10-S100 calcium binding protein A10 Gene

Reference: PTXBC015973

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	42C; ANX2L; ANX2LG; CAL1L; CLP11; Ca[1]; GP11; P11; p10; protein S100-A10; S100 calcium binding protein A10 (annexin II ligand, calpactin I, light polypeptide (p11)); annexin II ligand, calpactin I, light polypeptide; annexin II tetramer (AII) p11 subunit		

For research use only.

GENE DESCRIPTION	S100 calcium binding protein A10
GENE ID	6281
ORIGIN SPECIES	Human
ACCESSIONS	BC015973
VECTOR	pDONR223
NCBI SYMBOL	S100A10
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccatcctcaaatggaacacgccatggaaccatgatgtttacattcacaattcgtgggataaaggctacttaacaaaggaggacctgagagtcatggaaaaggagttccctggattttggaaaatcaaaaagaccctctggctgtggacaaaataatgaaggacctggaccag
tgtagagatggcgaatgggctccagagctcttttcctaatggggcctcaccattgcatgcaatgactattttgtagtacacatgaagcagaagggaagaagtag