

TOMM6-translocase of outer mitochondrial membrane 6 homolog (yeast) Gene

Reference: PTXBC015975

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	1110002E23Rik; AI663970; mitochondrial import receptor subunit TOM6 homolog; overexpressed breast tumor protein homolog; translocase of outer membrane 6 kDa subunit homolog; translocase of outer mitochondrial membrane 6 homolog (yeast)		

For research use only.

GENE DESCRIPTION	translocase of outer mitochondrial membrane 6 homolog (yeast)
GENE ID	100188893
ORIGIN SPECIES	Human
ACCESSIONS	BC015975
VECTOR	pDONR223
NCBI SYMBOL	TOMM6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcttcagcactgtcccgtgagcgtctggtcgtgtaataaactccgaaataccggacaacgtgggagattggttcgggctctaccgtttgccactgataaggaatgactccggaggaaactgatactaaattgggactcttgcgcgggagttggctggccaggaaac
ttgagtgcattgacctcatggcacctcagccaggggtgtag