



Reference: PTXBC015993

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein FAM151A; C1orf179; hyporthetical protein MGC27169; family with sequence similarity 151 member A		

For research use only.

GENE DESCRIPTION	family with sequence similarity 151, member A
GENE ID	338094
ORIGIN SPECIES	Human
ACCESSIONS	BC015993
VECTOR	pDONR223
NCBI SYMBOL	FAM151A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1