



Reference: PTXBC016029

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CTTNBP2 N-terminal-like protein; CTTNBP2 N-terminal like		
For research use only.			

GENE DESCRIPTION	CTTNBP2 N-terminal like
GENE ID	55917
ORIGIN SPECIES	Human
ACCESSIONS	BC016029
VECTOR	pDONR223
NCBI SYMBOL	CTTNBP2NL
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]