

CEP70-centrosomal protein 70kDa Gene

Reference: PTXBC016050

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BITE; centrosomal protein of 70 kDa; centrosomal protein 70kDa; p10-binding protein; testicular tissue protein Li 36; centrosomal protein 70		

For research use only.

GENE DESCRIPTION	centrosomal protein 70kDa
GENE ID	80321
ORIGIN SPECIES	Human
ACCESSIONS	BC016050
VECTOR	pDONR223
NCBI SYMBOL	CEP70
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtttccgtagccctaaacccagattccagtcacccatcagacagactcatgactgaaaaacagcaggaagaagcagaatgggaagcataaatgtgctattgatgcacatggcctaaacctttgtctctagtcaaaagaacagatctcaaagatctcatcttttgacaaacagtc
atcacaaaggatgagacagaatttgaaattgtgttggaagaacatcatgtcaacagaacatgatacaggagcttatagaactaatcaacagcttagaaatgaacttcagctagagcaaaagccgagcagccaatcaagaacaacgagctaatgacttggacaaattatggaaagtgtga
aatccaaattggtgaattggaggatgaatcactaaatagggcttgcaccaacagaataaaataaaagatctcaaaaggagcagaaaactttacaggtgaagtgcagcattataagaaaaaacgaacggagcaagaagaactattgtctcttgcacaaatggagctctagattaaaaa
aggaggaagaagatgcattgtcactcaaacagagtggttgcctatctgtcaaaagattcctcataccgtcttgatagacagtaa