

TMEM45B-transmembrane protein 45B Gene

Reference: PTXBC016153

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	transmembrane protein 45B		
For research use only.			

GENE DESCRIPTION	transmembrane protein 45B
GENE ID	120224
ORIGIN SPECIES	Human
ACCESSIONS	BC016153
VECTOR	pDONR223
NCBI SYMBOL	TMEM45B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcaattcaaggccacgcgctccagggagtttctctgacattgggctgtgttggtcagtgaaatgacccgctgaagtactttagccacacgcggaagaacagcccaactacattactatcagcgtctcagatcgtcgaagccgcaattaggactttgtttccgtcactgggatcctg
gcagagcagttgttccggatgggccccacctgcacctctaccatgagaaccactggataaagttaattgaattggcagcagacacacattgtacattcttgcagtcctcaggaaattgtgacatgctcacctatctggtcagccacgttccctgggggtggacagactgggtatggctgtggtc
agtattcatggaagggttctctctactaccacgtccacaacccgctcctggaccagcacatccactcactcctgctgtatgctctgttcggagggtgtgttagtatctccctagaggtgatcttccggaccacattgtgctggaactttccgaaccagctcactcattctcagggaaacct
ggttctggcagattgggtgtgtgctgttccaccttttggacacccgaatgggaccagaaggatgatgccaacctcatgttcacacatgtgcttctgctggcactacctggctgccctcagcattgtggccgtcaactattctctgttactgcttttgactcggatgaagagacacggaaag
gggagaatcatiggaattcagaagctgaattcagatgacattaccagaccgccctcttgagtggctcagatgaggaaatga