

LOH3CR2A-loss of heterozygosity, 3, chromosomal region 2, gene A Gene

Reference: PTXBC016278

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	LOH3CR2A; ERR-10; ERR10; NAG-7; NCRNA00312; long intergenic non-protein coding RNA 312		

For research use only.

GENE DESCRIPTION	loss of heterozygosity, 3, chromosomal region 2, gene A
GENE ID	29931
ORIGIN SPECIES	Human
ACCESSIONS	BC016278
VECTOR	pDONR223
NCBI SYMBOL	LOH3CR2A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcccatcattcttaaacacattttacatctggcacaataacgtgctacataccatctgtatttttctgccccatttttaaatcagccattctcaaggggtctcttcttaatctggctgtgtgtgtggaactcctgtaccatcttagaaccttaaggagacaagcgaaccaagccaataag
ctttccatgatgtgcttagagttaaacagagcccagggtactaagtattgtcatggagacagtgaaactaacctctggactgctgtctacatga