

BCL2A1-BCL2-related protein A1 Gene

Reference: PTXBC016281

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ACC-1; ACC-2; ACC1; ACC2; BCL2L5; BFL1; GRS; HBPA1; bcl-2-related protein A1; bcl-2-like protein 5; bcl2-L-5; hematopoietic BCL2-related protein A1; hemopoietic-specific early response protein; protein BFL-1; BCL2 related protein A1		

For research use only.

GENE DESCRIPTION	BCL2-related protein A1
GENE ID	597
ORIGIN SPECIES	Human
ACCESSIONS	BC016281
VECTOR	pDONR223
NCBI SYMBOL	BCL2A1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacagactgtgaatttgatatttacaggctggctcaggactatctgcagtgctctacagataccacaacctggatcagggtccaagcaaaacgtccagagtgtacaaaatgttgcgttctcagtcacaaaagaagtggaaaagaatctgaagtcattgtggacaatgttaattgtgt
gtccgttagacactgccagaacactattcaaccaagtgtggaaggagttgaaagacggcatcatttaactggggaagaattgtaacctatttgcatttgaaggtattctcatcaagaacttctacagacaaattgccccggatgtggatacctataaggagatttcatatttgttgcggag
ttcataatgaataacacaggagaatggataaggcaaacggagctgggaaatggctttgtaagaagtgtgaacctaaatctggctggatgacttttctagaagttacaggaaagatctgtgaaatgctatctctctgaagcaatactgttga