

CXCL3-chemokine (C-X-C motif) ligand 3 Gene

Reference: PTXBC016308

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CINC-2b; GRO3; GROg; MIP-2b; MIP2B; SCYB3; C-X-C motif chemokine 3; GRO-gamma; GRO-gamma(1-73); GRO3 oncogene; MGSA gamma; MIP2-beta; chemokine (C-X-C motif) ligand 3; growth-regulated protein gamma; macrophage inflammatory protein 2-beta; melanoma growth		

For research use only.

GENE DESCRIPTION	chemokine (C-X-C motif) ligand 3
GENE ID	2921
ORIGIN SPECIES	Human
ACCESSIONS	BC016308
VECTOR	pDONR223
NCBI SYMBOL	CXCL3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcccaagccacgtctctccgccccccagcaatccccggctcctgcgggtggcgtgctgctcctctctctggtggccgccagccggcgccagcaggagcgtccgtggtcactgaactgcgtgccagtgtcttcagacactgcagggaattcacctcaagaacatccaaagtgtgaatgtaaggtcccccgaccacctgcgccccaaaccgaagtcatagccacactcaagaatgggaagaagctgtctcaaccccgcatccccatggttcagaaaatcatgaaaagatactgaacaaggaggagcaccaactga