

C1orf149-chromosome 1 open reading frame 149 Gene

Reference: PTXBC016328

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C1orf149; CENP-28; EAF6; NY-SAR-91; chromatin modification-related protein MEAF6; Esa1p-associated factor 6 homolog; centromere protein 28; protein EAF6 homolog; sarcoma antigen NY-SAR-91; MYST/Esa1 associated factor 6		

For research use only.

GENE DESCRIPTION	chromosome 1 open reading frame 149
GENE ID	64769
ORIGIN SPECIES	Human
ACCESSIONS	BC016328
VECTOR	pDONR223
NCBI SYMBOL	C1orf149
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgtgcacaacaaggcggcgccgcgcagatcccggacacccggcgggagctggcggagctcgtgaagcggagcaggagctggcggaaacattggcaatttgagcgcagatctatgcttttggaggagctacgtggaagacactcagatgtatggcaatatttcgt
ggctgggatcggtatctgaccaaccccccaactccaatagcaaaaatgatcgaaggaaacgggaagttaagggaagctgagcggctctcagtaaatcctcggttacctcagcagctgcagtaagtgcattggcaggagttcaggaccagctcattgaaaaggaggcagggaagtgg
gacggaaagtgcacttctccagacttccacaatcaggaaaatgagccagccaggagaccctgaggatctggtggtatctgtgcaggagtgaaacctcagaaggctgcttcttctacttctcagggagtcaccacagcagccataaaaagcgaagaataaaaccggcacagg
attgatctgaagttaacaaaaaacacagctgactattag