



Reference: PTXBC016440

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	IOP2; nuclear prelamin A recognition factor; iron-only hydrogenase-like protein 2; prenyl-dependent prelamin A binding protein		

For research use only.

GENE DESCRIPTION	nuclear prelamin A recognition factor
GENE ID	26502
ORIGIN SPECIES	Human
ACCESSIONS	BC016440
VECTOR	pDONR223
NCBI SYMBOL	NARF
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]