

EDEM3-ER degradation enhancer, mannosidase alpha-like 3 Gene

Reference: PTXBC016464

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	alpha-1,2-mannosidase EDEM3; C1orf22; ER degradation-enhancing alpha-mannosidase-like protein 3; ER degradation enhancer, mannosidase alpha-like 3; ER degradation-enhancing -mannosidase-like protein 3; ER degradation-enhancing alpha-mannosidase-like 3; ER		

For research use only.

GENE DESCRIPTION	ER degradation enhancer, mannosidase alpha-like 3
GENE ID	80267
ORIGIN SPECIES	Human
ACCESSIONS	BC016464
VECTOR	pDONR223
NCBI SYMBOL	EDEM3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggaggcttaagtatggtgatgaacaaccacgttctgactggcgtagtatcgtaggaatctggagcatgctgtgttagaattgacctgtttaaactgtccatcaaaaatggaaatccacagttccccctcaaatgcagcactgcaccacctgcaacacctcaggccagggaaagat
tactgagcattctcgcgaaccagatttctgttctctgtagacaagaacaataatcatttagtagtgagtgagggaataggagtttgtagccttctaattaaaggagctgccttcttggtttgggggttagaggctctttgggaagatgcacccgagtagttgtggcattctga