

## HHIPL2-HHIP-like 2 Gene

Reference: PTXBC016552

|                               |                                                                                  |                     |                |
|-------------------------------|----------------------------------------------------------------------------------|---------------------|----------------|
| Product type:                 | DNA                                                                              | Storage Condition:  | -20°C          |
| Size/Quantity:                | 2 ug                                                                             | Delivery Condition: | +4°C, Blue Ice |
| Synonyms                      | KIAA1822L; HHIP-like protein 2; hedgehog interacting protein-like 2; HHIP like 2 |                     |                |
| <b>For research use only.</b> |                                                                                  |                     |                |

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | HHIP-like 2                                                                                      |
| <b>GENE ID</b>          | 79802                                                                                            |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC016552                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | HHIPL2                                                                                           |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon)

atggtggcccaagatcctcagggtgctgcagctgcctgagcgaggtggccacgggctgaggaaacccgctcctcatggtccatgctggggacggcaccatcctcttctgttggccgagcaggtaggagtggtgggtctacctccctgatgggagtcgctggagcaacccttc  
ctggacctcaagaacatcgtgttgaccacccatggatcggggatgagagaggtcttgggttggcttttaccaccaatccgccacaatcgcaagtctctatttatttctgctggacaagaagaaggtagaaaagatccgaattagagatgaaggtttctcgggtgatcctaa  
caagctgacctgaaatcagagagggtcatcttggaattgaagaaccagctcaaacataatggcggacaactcttttggcctggtatgctatatgtacatattcactggggacgggggacagctggagatcccttggcctgttggaaatgctcagaacaaagtccctgctggg  
aaaagtgttaaggatcgtggaacaggggcaggctcacatggcaagcgttacggagtcctcctgggacaatccattgtttctgagccaggggccaccgccatctatgcctatgggatcaggaaacatgtggcgttctgtgtgacgggggaccccatcacgcgcaggggccgag  
gccgatatctgtggggacgtggggccagaacaggttgaagaggttgacctatgttgaagggtgaaactatggctggagagcaaaaggaggttgcattgtatgacaaaaactttgcacatgctcttggatgatgtctgccaatctatgcttatggccatgagtgaggagca  
gtcactggaggttatgtctatcgtgttgaatcccaaatcgaatggcctgtatcttggagactcatgagtgctgacttatgcttgcaggaaagatagaaaaacaagaatggaagaagcaggatcttgcctgggcagcaccacgtcctgtgcttccagggtgatcagcac  
ccatagcaagttcatcatctcttctgaagatgaagcaggggagctgtatcttggcgacctcttaccgaagtgcctatgcaccagtgatctattacaagtgttggacccctcaaggcgagcaccgcccaaggcaagtgcaaatcaagccagtgccgtgagaaccaagagtaagcg  
gatccgttcagaccactgcgaagacagcttggacttgaaggaacaatcagaagaagctgtgagaaatctccagtgcacaccttagcttctgcccagcccagggttctgtgagaaaggctcctccaagaagctggtctctcacaagcagcaagaatactgcagggcctg  
gtacaaagaagaagccagagtggtggggcccccacgtccgccaggcgaagaggaggaagagcctgaaaagccacagtggtgcaggatgagccatcagcagcagaagcgagctggcagaagctccttga