



Reference: PTXBC016663

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	VDU1; ubiquitin carboxyl-terminal hydrolase 33; VHL-interacting deubiquitinating enzyme 1; deubiquitinating enzyme 33; pVHL-interacting deubiquitinating enzyme 1; ubiquitin thioesterase 33; ubiquitin thiolesterase 33; ubiquitin-specific-processing proteas		

For research use only.

GENE DESCRIPTION	ubiquitin specific peptidase 33
GENE ID	23032
ORIGIN SPECIES	Human
ACCESSIONS	BC016663
VECTOR	pDONR223
NCBI SYMBOL	USP33
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]