

THNSL1-threonine synthase-like 1 (*S. cerevisiae*) Gene

Reference: PTXBC016697

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TSH1; threonine synthase-like 1; threonine synthase like 1		
For research use only.			

GENE DESCRIPTION	threonine synthase-like 1 (<i>S. cerevisiae</i>)
GENE ID	79896
ORIGIN SPECIES	Human
ACCESSIONS	BC016697
VECTOR	pDONR223
NCBI SYMBOL	THNSL1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatgggaatcccgattcgaaaatttatctgtgcctctaatacagaaccatgttttgactgattttataaaaacaggacattatgatctaagggaagaaaactagcacaaacctttaccgtcaatagatatctcaaatctcaaacctagaacgacatttacacttgatggcctaataagatggac
agctaafgacagaattatttaatacattagaaaagtcagcatctccagatagaaaaggtctagttgagaaacttcagcaggtattttgtagctgactgggtgctctgaggagagtgccctagcagctattaactccactataatactcagggatatttggatccacactgctgttgcataaa
gtggttgcatagagggtgcaagacaaaactgcccgtgattatctcatctacagcccactactcaagtttgcacctgctatcatgcaggctttaagattaaagaaatcaatgagactcatcaagtcagctctatttctgggttcatacaatgcattacccactgcattgaggctttattagaga
gaacaaaacagcaagagaagatggagtaccaggtctgtgcagctgatatgaatgtcttgaagagtcgtggaacaactgtccaaaatcaattcatatga