

CHMP5-chromatin modifying protein 5 Gene

Reference: PTXBC016698

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C9orf83; CGI-34; HSPC177; PNAS-2; SNF7DC2; Vps60; charged multivesicular body protein 5; SNF7 domain containing 2; SNF7 domain-containing protein 2; apoptosis-related protein PNAS-2; chromatin-modifying protein 5; hVps60; vacuolar protein sorting-associated		

For research use only.

GENE DESCRIPTION	chromatin modifying protein 5
GENE ID	51510
ORIGIN SPECIES	Human
ACCESSIONS	BC016698
VECTOR	pDONR223
NCBI SYMBOL	CHMP5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaccgactctcgggaaagcgaaccacaaaggctccgccgccagcctgactggctgcaatggcaggtggacagtagagcagaatccattgacaagaagattctcgattggatgctgagctagtgaagtataaggatcagatcaagaagatgagagagggctcgtcaagaatat
ggfcaagcagaagccttcgaggttttaagcaaaagaggatgatgagcagcagcgggacaattctgcccaacagtcattcaacatggaacaagccaattatccatccagctcttgaaggacaccaagaccagcgttggatgctatgaaactgggagtaaaagaaatgaagaaggcata
caagcaagtgaagatcgaccagattgaggattacaagaccagctagaggatgatggaagatgcaaatgaaatccaagaagcactgagtcgagttatggcaccacagaactggatgaagatgattagaagcagagttggatgcactaggtgatgagcttctggctgatgaagacagt
tcttatttgatgaggcagcatctgcacctgcaattccagaaggtgttcccactgatacaaaaaacaaggatggagttctgtggatgaattggattgccacagatccctgctcatag