

SH3BGR1-SH3 domain binding glutamic acid-rich protein like Gene

Reference: PTXBC016709

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HEL-S-115; SH3BGR; SH3 domain-binding glutamic acid-rich-like protein; SH3 domain binding glutamic acid-rich protein like; SH3-binding domain glutamic acid-rich protein like; epididymis secretory protein Li 115; SH3 domain binding glutamate rich protein 1		

For research use only.

GENE DESCRIPTION	SH3 domain binding glutamic acid-rich protein like
GENE ID	6451
ORIGIN SPECIES	Human
ACCESSIONS	BC016709
VECTOR	pDONR223
NCBI SYMBOL	SH3BGR1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtgatccgtgtatatattgcattctctctggtctacagcgattaaagaagaacaagaatgtgcttggttcctagaagccaacaaataggattgaagaaaagatattgcagccaatgaagagaatcggaaaggatgagagaaaatgtacctgaaaatagtcgaccagccacaggttaccctgccacctcagattttcaatgaaagccagatcgcgggactatgatgcctcttgaagccagagaaaaatagcagtgatgcctcttaggcttgacagccccacctggttcaagggaagcagaagtgcagcaagcaagcagcaagcatga