

## HAX1-HCLS1 associated protein X-1 Gene

Reference: PTXBC016730

|                |                                                                                                                                                                                                                  |                     |                |
|----------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                                              | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                                             | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | HCLSBP1; HS1BP1; SCN3; HCLS1-associated protein X-1; HAX-1; HCLS1 (and PKD2) associated protein; HS1 binding protein; HS1-associating protein X-1; HS1-binding protein 1; HSP1BP-1; HCLS1 associated protein X-1 |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | HCLS1 associated protein X-1                                                                      |
| <b>GENE ID</b>          | 10456                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC016730                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | HAX1                                                                                              |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atgagcctcttgcctcctcggggcttttgccttctgacctcggagccacagagatccctttttggaggggactcagatgaagatgatgatgaggaagaagaagaagggggctcatggggccgtgggaaccaaggttcatagtcctcagacccccctgaggaattt  
ggcttcggcttcagcttcagccaggaggaggatagctttccacgataacttcggctttgatgacgtacgagattcaatagcatcttcagcgatatggggcctggaccttgcttccatcctcctgaactccaggtcctgagtcggagacacctggtgagagactacggaggga  
cagacacttcgggactcaatgcttaagtatccagatagtcaccagccaggatctttggggggctcttgagagtgatgcaagaggtgaatcccccaaccagcaccagactggggctccagaggccattcataggtttgatgatgatggcctatggaccccatcctagaaccagag  
aggacaatgatcttgattccaggtttccaggagggtcttggccgggttctacagccagcccaatcctattcaagagcatctctgtgaccaagatcactaaaccagatgggatagtggaaggagcggcggactgtgtgacagtggggccggacagagactacagtaaccgac  
acgaagcagatagcagtcctagggtgatccagaatcaccaagacctccagccctggatgatgcttttccatcctggacttattcctgggacgttgggttcgggtcccgtag