

TMSB10-thymosin beta 10 Gene

Reference: PTXBC016731

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MIG12; TB10; thymosin beta-10; migration-inducing gene 12; migration-inducing protein 12; thymosin beta 10		

For research use only.

GENE DESCRIPTION	thymosin beta 10
GENE ID	9168
ORIGIN SPECIES	Human
ACCESSIONS	BC016731
VECTOR	pDONR223
NCBI SYMBOL	TMSB10
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcagacaaccagacatggggaaatcgccagctcgataaggccaagctgaagaaaacggagacgcaggagaagaacacccctgccaccaagagaccattgagcaggagaagcggagtgaaatttctaa