

**BRF1-BRF1 homolog, subunit of RNA polymerase III transcription initiation factor IIIB
(S. cerevisiae) Gene**
Reference: PTXBC016743

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BRF1, RNA polymerase III transcription initiation factor subunit; BRF1, RNA polymerase III transcription initiation factor 90 kDa subunit; BRF1 homolog, subunit of RNA polymerase III transcription initiation factor IIIB; BRF; BRF-1; CFDS; GTF3B; HEL-S-76p		

For research use only.

GENE DESCRIPTION	BRF1 homolog, subunit of RNA polymerase III transcription initiation factor IIIB (S. cerevisiae)
GENE ID	2972
ORIGIN SPECIES	Human
ACCESSIONS	BC016743
VECTOR	pDONR223
NCBI SYMBOL	BRF1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacgggcccgcgtgtccgcggttgcggcgacggacatcgagctggacggcgccggggacggcgtgtgcaccgcctgcggcagtcgtggaggacaacatcgtgtccgaggtgcagttcgtggagagcagcggcgccgctcctcggccgtggccagttcgtgtccctggacggtgtgtgcaaaacccgactctgggtggcggttccacgtgaatctggggaaggagtcgagagcgcagaccctgcagaatggagggcggccacatccaccacctggggaaccagctgcagctgaaccagcactgcctggacaccgcttcaacttcaagatggccgtgagcaggcaccctgacccggccgggaagatggccacgtgattgtgcctgcctctacgtgtgtcgtacggaggcagcggcgcacatgtcctggacctcagcgacctgctccaggtagacagcctccgtcctgcacttccccacctggggtgtgacctgggggtgtgaccagggtgtgaccgggtgtaccacaggtgcctccacgcacatcagtgcccggtctgtgtgcctgcccggtcagggaagtgtgtgtgttagatga