

C4orf36-chromosome 4 open reading frame 36 Gene

Reference: PTXBC016746

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C4orf36; chromosome 4 open reading frame 36		
For research use only.			

GENE DESCRIPTION	chromosome 4 open reading frame 36
GENE ID	132989
ORIGIN SPECIES	Human
ACCESSIONS	BC016746
VECTOR	pDONR223
NCBI SYMBOL	C4orf36
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgtatggagtgccaagaagaacacagtgaaaaccatttgcggggcagttgtataatgtacaggaaccttgggatatgcatgctgcaagacctgggtacacaaacctagccaatatcaagttgccttcttggagaatttcatttgggtgtctgtgcagctcacaaatgtacca
ccattaaagatggactgtccctcttcagaaatctcaaaactgaaaggagatgaagtgaagcgtcttgaactgaagtgtcaagaaaatacatctaaggaattcagcttctcctgagggaaaggccagccggtttgagaagacctcttccatctaaatga