



Antibodies and Proteins

## ProteoGenix SAS

## Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

[www.ProteoGenix.science](http://www.ProteoGenix.science)

## PDLIM4-PDZ and LIM domain 4 Gene

Reference: PTXBC016765

|                |                                                                                                                                        |                     |                |
|----------------|----------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                    | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                   | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | PDZ and LIM domain protein 4; LIM domain protein; LIM protein RIL; enigma homolog; reversion-induced LIM protein; PDZ and LIM domain 4 |                     |                |

**For research use only.**

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | PDZ and LIM domain 4                                                                              |
| <b>GENE ID</b>          | 8572                                                                                              |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC016765                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | PDLIM4                                                                                            |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

**ORF (Only for reference, product provided without STOP codon )**

atgcccaatcctgtagaccctgcgcgggcttcgccctcggggcttcgccctgctgctggggcgccgggactcagcgcgccctaccatcaccagggttcacatcgtgcgagcaaggctcgaattggctgcctctgcccaggagacctgatccaggccatcaatggtgagagcacagagctcat  
gacacacccctggggcgacagcaagccgaatcaaggctgcgcacgatcactcactacatctctgtgagcagacgtcaggcgaggagctcgcgcggccttcagctgatgacagacagcctcagcagacacagatcagatcagatcagacgagcgagcccaacacacagagctg  
gccctcaccggcgccagactggcccgagaagatgtagcagcaaacgctggatctctatgtagcaaacccctcgtcttccagtcctccatcgaatggggcgccgcctcagggccaccctcagggccagatgacacctgatgctctccaccccagcctgaccacccagagagcctccc  
cgcgagcgggagctcagagctcgaactcgggtctcgagggtgtacagagctgctcgggagcagccagagccctgctggccgcggagcccaagcagtcagcgtctctccgctactctgagggcatctagagggccgcggagccggggattggcccgcccgccgggac  
ctcaagccacccagcgcagcaagctcctcctctcagcgcctcagggggctcccgagctgcagcgcctgcgcggccacgcatctggcgccacatctgagccagcagcgggacagcctaccatcccgagtgctctcatctgctcagtgactcggcctgaaactcaagcagcgtgttta  
ctcttctctgcagcgaagcgtctactctctgtagcagcgcgcgaagcgcgcgtctgaagcccccagagctcagcagctgctcggctgtatccccaatgcgaagctggaactcctctga

More info: Please contact us for any question: [contact@proteogenix.fr](mailto:contact@proteogenix.fr)