



Reference: PTXBC016852

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HP1P; pre-B-cell leukemia transcription factor-interacting protein 1; hematopoietic PBX-interacting protein; pre-B-cell leukemia homeobox interacting protein 1; PBX homeobox interacting protein 1		

For research use only.

GENE DESCRIPTION	pre-B-cell leukemia homeobox interacting protein 1
GENE ID	57326
ORIGIN SPECIES	Human
ACCESSIONS	BC016852
VECTOR	pDONR223
NCBI SYMBOL	PBXIP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr