



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

ATF4-activating transcription factor 4 (tax-responsive enhancer element B67) Gene

Reference: PTXBC016855

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CREB-2; CREB2; TAXREB67; TXREB; cyclic AMP-dependent transcription factor ATF-4; DNA-binding protein TAXREB67; cAMP response element-binding protein 2; cAMP-dependent transcription factor ATF-4; cAMP-responsive element-binding protein 2; cyclic AMP-respon		

For research use only.

GENE DESCRIPTION	activating transcription factor 4 (tax-responsive enhancer element B67)
GENE ID	468
ORIGIN SPECIES	Human
ACCESSIONS	BC016855
VECTOR	pDONR223
NCBI SYMBOL	ATF4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr